



Environmental Fast BTEX Analyzers



www.dps-instruments.com

When you need speedy BTEX (Benzene, Toluene, Ethylbenzene, and Xylenes) analyses to satisfy your clients choose one of our full range of BTEX Analyzers. We offer both Series 600 Lab GC's and Companion Portable GC's configured for Direct Injection for your extracts, or gas samples. For the automated analysis of gas samples, choose the Gas Sample Valve option. The integrated Sample loop is either filled using a sample bag, or an electronically controlled vacuum pump can automatically draw the sample through the Sample Loop for highly precise and reproducible analyses. All of the BTEX Analyzers are fully integrated GC Systems that are small and lightweight and all DPS systems are modular for expandability, upgrades, and easy service.




Series 600 GC

Available Configurations Include:

- 600-C-010 - Series 600 BTEX GC Analyzer (FID, PID, 30m)
- 500-C-010 - Companion 1 Portable BTEX GC Analyzer (PID, 30m)

100 ng BTEX with PID and FID in Series



11/2019
Specifications may change without notice.

Companion 1 Portable GC
(with Gas Sample Valve)